



RIDA®CUBE SCAN

The smart analyzer

For detailed
information
please
follow us



RIDA® CUBE SCAN

Three analysers are available with different combinations of wavelengths: 340 nm is always in place, in combination with a second wavelength respectively at 505, 546 or 580 nm. The combination of different instruments allows to increase the throughput and to cover all wavelengths for the future test-kit portfolio. The analyser set includes the instrument and the tablet PC. The tablet PC is also available separately if broken or lost.

Technical data

RFID technology

- Pre-set calibration curve
- Lot number and expiry
- Number of tests remaining

Maintenance free

- Fully-automated system without pipetting device
- Closed test-cartridges
- Integrated QC test

Tablet PC interface

- Android based app, updates available via Internet
- Bluetooth/USB connection
- Quick data entry with pictograms
- Feed-back option for trouble shooting
- Data transfer to host printer

Specifications

- Weight: 2.4 kg
- Dimensions: 16 x 13 x 14.5 cm
- CE-compliant